








| Item  | Features   | Procurement   |
|---|--|---|
| <p><b>SoftEdge™</b></p>                            | <ul style="list-style-type: none"> <li>Made with an expanded-vinyl foam that's extra-soft, reducing pressure on skin while protecting you from hard work-surface edges. The softest edge protectors available. Choose yellow to alert and protect workers. Can be cut into custom lengths. Smooth surface wipes clean. Hook and loop fastener for secure attachment.</li> <li>2½"W x ½" thick, 1¾" side</li> </ul> | <p><a href="#">AliMed</a></p>                         |
| <p><b>SoftEdge™<br/>Corners</b></p>               | <ul style="list-style-type: none"> <li>Covers sharp corners and make continuous Softedge protection possible (4 pack)</li> </ul>   | <p><a href="#">AliMed</a><br/>#75189<br/>(4 pack)</p> |
| <p><b>AliEdge™ Gel Edge<br/>Protector 1</b></p>  | <ul style="list-style-type: none"> <li>Provides a soft, completely sealed buffer between employees and hard, sharp edges</li> </ul>  | <p><a href="#">AliMed</a><br/>#792307</p>             |
| <p><b>AliMed®<br/>Ulnar Gel Pads</b></p>         | <ul style="list-style-type: none"> <li>Gel pads to place under elbow</li> <li>Sold in pairs</li> <li>Relieve contact pressure on the ulnar nerve</li> <li>Useful for pipetting and other high-precision tasks where the elbow rests on a hard surface</li> </ul>   | <p><a href="#">AliMed</a><br/>#78336</p>              |

|  |   |   |
|--|---|---|
| <p><b>Wedge-Ease Ergonomic Forearm Supports</b></p>     | <ul style="list-style-type: none"> <li>• Vinyl w/ nonslip base</li> <li>• Four different sizes can be combined/stacked</li> <li>• Also available in ESD</li> <li>• Rated Fed Std. 209 and ISO14644 Clean Room Class 100 III</li> <li>• Can be used to provide both flat and angled support</li> <li>• Helpful for use with microscope work to reduce neck and shoulder muscle fatigue</li> <li>• Reduces the contact pressure at forearms and wrists</li> <li>• Item #79597 3" wedge</li> </ul> | <p><a href="#">AliMed</a><br/>Item #79598<br/>2" Wedge</p> <p>Item #79599<br/>Regular<br/>Stacker</p> <p>Item #79600<br/>Mini Stacker</p> |
| <p><b>Morency Forearm Support for Microscopes</b></p>  | <ul style="list-style-type: none"> <li>• 5" and 8" wide versions</li> <li>• Includes a clamp that fits 5/8" to 1 3/4" thick work surfaces</li> <li>• Adaptor available for thicker workstations</li> <li>• Angle-adjustable forearm support for microscope users</li> </ul>   | <p><a href="#">Morency</a><br/>#MSM-N-BL-T<br/>5" wide</p> <p>#MSM-N-BL<br/>8" wide</p>   |
| <p><b>Deluxe Edge Rest</b></p>                        | <ul style="list-style-type: none"> <li>• Attaches securely with self-adhesive hook and loop tabs.</li> <li>• 1 1/8" thick expanded urethane</li> <li>• Black. 22" L x 5" W</li> <li>• Self-skinned foam surface wipes clean.</li> <li>• Fits over top edge of table surface for protection from sharp edges</li> </ul>  | <p><a href="#">AliMed</a><br/>Item #73075</p>   |
| <p><b>Biohit</b></p>                                  | <ul style="list-style-type: none"> <li>• Peel and stick pads can be cut down to size</li> <li>• Closed cell foam SoftKnees pads can be cut for elbows &amp; forearm pads to fit inside sleeves</li> <li>• Universal Peel and Stick Pads</li> </ul>  | <p><a href="#">Softknees</a><br/>#1012R</p>   |

**Nuaire Foam Armrest pads**



- Disposable white foam pad
- 3/16" closed cell foam
- Does not absorb liquid
- Package of 10
- Self-stick adhesive tape
- Easily removed and replaced
- Easy to remove and replace foam elbow pads used to pad surface of biological safety cabinet to reduce soft tissue compression and discomfort

[Scientific  
Visions](#)  
#NU-976